

EMIT ANIMAL AND TISSUE PREPARATION PROCEDURE

Pre-preparation Setup

Materials

Whole Mouse (Preferred Method)	Whole Mouse (Method 2) and Tissues
Dewar flask	Cooler and Metal pan
Hexanes, ACS grade	Aluminum foil
Dry ice	Dry ice
Clear plastic bags	Clear plastic bags
Sharpie/Permanent marker	Sharpie/Permanent marker



Dewar



Cooler



Pan

- Euthanize animals using CO₂.
- Begin the freezing procedure immediately following euthanasia.
- It is recommended that all steps are performed in the same room

Freezing Procedures

I. Freezing whole mice – Preferred Method

Note: the preferred method for freezing whole animals is submersion in a cold bath of hexane and dry ice (isopentane may be substituted for hexanes). If the use of hexanes is not permitted in your facility, use Method 2. We **do not** recommend freezing animals in liquid nitrogen as it freezes too quickly and can result in cracks.

Note: Hexanes are volatile and flammable, perform the steps in a fume hood, and wear safety glasses and thermal gloves.

1. Place dry ice in a tall Dewar flask, using enough to cover the bottom.
2. Slowly add the hexanes to the dewar.
 - Add hexanes to a height that will allow full submersion of the mouse without touching it to the dry ice at the dewar bottom.
3. Allow the temperature of bath to equilibrate to the dry ice freezing point (-78°C).
 - Boiling of the hexanes will cease once the bath temperature is equilibrated.
 - There should still be dry ice at the bottom of the dewar, add more if needed.
4. Hold the mouse by the tail and slowly lower it into the cold bath.
 - The hexanes should not boil over during this process.
 - It should take a few minutes to submerge the entire mouse.
5. Once the mouse is fully submerged, let it remain in the hexanes for 10 additional minutes.
6. Remove the mouse from the bath and store individually in a clearly labeled bag.
7. Place the sample in the freezer and store at -80°C.

II. Freezing whole animals – Method 2:

1. Cover the bottom of a cooler with crushed dry ice.

2. Place an aluminum foil-lined metal pan on the crushed dry ice.
3. Cover the cooler and allow at least 10 minutes for temperature equilibration.
4. Place the mouse on the foil-lined pan, belly side down, with arms and legs extended.
5. Quickly cover the mouse with **powdered** dry ice.
6. Cover the cooler and leave for 45 minutes to allow the mouse to freeze solid.
7. Once frozen, transfer the mouse to a clearly labeled bag and store at -80°C.

III. Freezing tissue specimens

1. Cover the bottom of a cooler with crushed dry ice.
2. Place an aluminum foil-lined metal pan on the crushed dry ice.
3. Cover the cooler and allow at least 10 minutes for temperature equilibration.
4. Place the tissue in the desired conformation/orientation onto the foil-lined pan.
5. Cover the cooler and leave for 5-10 minutes to allow the tissue to freeze in its natural state.
6. Wrap the tissue in the foil that it was chilled on.
 - **Note:** *Tissue must be frozen prior to wrapping! Wrapping up the tissue before it is frozen will destroy the anatomical morphology.*
7. Wrap all samples carefully in additional pre-cooled, labeled foil.
8. Place in a pre-chilled, labeled plastic bags.
9. Place the samples in a pre-chilled freezer box and store at -80°C.

If you need any further assistance or have questions, contact info@emitimaging.com